Laboratory	Central Forensic Science Laboratory, 30, Go West Bengal	rachand Road	l, Kolkata,
Accreditation Standard	ISO/IEC 17025: 2005		
Discipline	Forensic Testing	Issue Date	20.11.2015
Certificate Number	T-1576	Valid Until	19.11.2017
Last Amended on	-	Page	1 of 17

S. No.	Product / Material of Test	Specific Test Performed	Test Method Specification against which tests are performed	Range of Testing / Limits of Detection
I.	BIOLOGY & DNA			
1.	Blood	Tetramethyl Benzidine Test	Working Procedure Manual in Biology; CFSL/KOL/BIO/WPM Issue No. 01; Date of Issue: 01.01.2007	106 dilution
		Takayama Test	Working Procedure Manual in Biology;CFSL/KOL/BIO/WPM Issue No. 01; Date of Issue: 01.01.2007	1 μl of Blood or 0.1 mg of Haemoglobin
		Anti Human Haemoglobin Test	Working Procedure Manual in DNA Profiling: CFSL/KOL/DNA/WPM Issue No. 01; Date of Issue: 01.01.2007	0.05 μg of Haemoglobin
2.	Semen	Acid Phosphatase Test	Working Procedure Manual in Biology;CFSL/KOL/BIO/WPM Issue No. 01; Date of Issue: 01.01.2007	1 μl of Semen
		P30 Test	Working Procedure Manual in DNA Profiling: CFSL/KOL/DNA/WPM; Issue No. 01; Date of Issue: 01.01.2007	1 μl of Semen
		Microscopic Examination	Working Procedure Manual in Biology; CFSL/KOL/BIO/WPM Issue No. 01; Date of Issue: 01.01.2007	0.5 µl of Semen

Laboratory	Central Forensic Science Laboratory, 30, Gorachand Road, Kolkata, West Bengal		
Accreditation Standard	ISO/IEC 17025: 2005		
Discipline	Forensic Testing	Issue Date	20.11.2015
Certificate Number	T-1576	Valid Until	19.11.2017
Last Amended on	-	Page	2 of 17

S. No.	Product / Material of Test	Specific Test Performed	Test Method Specification against which tests are performed	Range of Testing / Limits of Detection
3.	Saliva	Starch Iodine Test	Working Procedure Manual in Biology; CFSL/KOL/BIO/WPM Issue No. 01; Date of Issue: 01.01.2007	10 μl of Saliva
		Saligae Test	Working Procedure Manual in DNA Profiling: CFSL/KOL/DNA/WPM Issue No. 01; Date of Issue: 01.01.2007	5 μl of Saliva
4.	Hair	Morphologic and Microscopic Examination	Working Procedure Manual in Biology; CFSL/KOL/BIO/WPM Issue No. 01; Date of Issue: 01.01.2007	1 strand for origin; >10 strands for comparison
5.	Blood for species identification	Double Diffusion Test	Working Procedure Manual in Biology; CFSL/KOL/BIO/WPM; Issue No. 01; Date of Issue: 01.01.2007	> 1µL of Blood
6.	Human Blood, Semen, Vaginal secretions, Saliva, Bone, Tooth, Muscle, Viscera & other	DNA Extraction	Working Procedure Manual in DNA Profiling: CFSL/KOL/DNA/WPM Issue No. 01; Date of Issue: 01.01.2007	>500 nucleated cells
	Tissues	Yield Gel	Working Procedure Manual in DNA Profiling: CFSL/KOL/DNA/WPM Issue No. 01;Date of Issue: 01.01.2007	20 ng to 100 ng

Laboratory	Central Forensic Science Laboratory, 30, Gorachand Road, Kolkata, West Bengal		
Accreditation Standard	ISO/IEC 17025: 2005		
Discipline	Forensic Testing	Issue Date	20.11.2015
Certificate Number	T-1576	Valid Until	19.11.2017
Last Amended on	-	Page	3 of 17

S. No.	Product / Material of Test	Specific Test Performed	Test Method Specification against which tests are performed	Range of Testing / Limits of Detection
	Human Blood, Semen, Vaginal secretions, Saliva, Bone, Tooth, Muscle, Viscera & other	Amplification of Short Tandem Repeats (Autosomal & Y Chromosomal)	Working Procedure Manual in DNA Profiling: CFSL/KOL/DNA/WPM Issue No. 01; Date of Issue: 01.01.2007	0.5 ng to 1 ng
	Tissues	Fragment Sizing and Genotyping with software	Working Procedure Manual in DNA Profiling: CFSL/KOL/DNA/WPM Issue No. 01; Date of Issue: 01.01.2007	(±)0.4 to 0.5 bases
П.	BALLISTICS	Statistical Analyses	Working Procedure Manual in DNA Profiling: CFSL/KOL/DNA/WPM Issue No. 01; Date of Issue: 01.01.2007	Qualitative
1.	Firearms	Physical Examination Test Firing Chemical Analysis	Working Procedure Manual In Ballistics; CFSL/KOL/BALL/WPM Issue No. 01; Date of Issue: 01.01.2007	Qualitative
2.	Pellets, Bullets, Cartridge Cases, Cartridges and its Components	Physical Examination Microscopic Examination Test Firing Comparative Examination	Working Procedure Manual In Ballistics; CFSL/KOL/BALL/WPM Issue No. 01; Date of Issue: 01.01.2007	Qualitative

Laboratory	Central Forensic Science Laboratory, 30, Gorachand Road, Kolkata, West Bengal		
Accreditation Standard	ISO/IEC 17025: 2005		
Discipline	Forensic Testing	Issue Date	20.11.2015
Certificate Number	T-1576	Valid Until	19.11.2017
Last Amended on	-	Page	4 of 17

S. No.	Product / Material of Test	Specific Test Performed	Test Method Specification against which tests are performed	Range of Testing / Limits of Detection
3.	Targets of Firing	Physical Examination	Working Procedure Manual In Ballistics;	
		Microscopic Examination	CFSL/KOL/BALL/WPM Issue No. 01;	Qualitative
		Test Firing	Date of Issue: 01.01.2007	
		Comparative Examination		
		Instrumental Examination		
4.	Gun Shot Residue	Physical Examination	Working Procedure Manual In Ballistics;	Qualitative
		Chemical Analysis	CFSL/KOL/BALL/WPM Issue No. 01;	
		SEM-EDX Analysis	Date of Issue: 01.01.2007	
III.	FORENSIC ELEC	TRONIC INVESTIGATION (CO	OMPUTER FORENSIC)	
1.	Hard disk	Disk forensic ; Bit stream imaging, cloning	Working Procedure Manual in Computer Forensic CFSL/KOL/CF/WPM	Qualitative
		OS analysis, Data analysis	Issue No. 01; Date of Issue: 01.01.2007	
		Internet artifacts, Registry forensic, E-mail Social		
		networking site artifacts and chat extraction.		
		Retrieval of detected / hidden data / timeline examination and keyword search		

Laboratory	Central Forensic Science Laboratory, 30, Gorachand Road, Kolkata, West Bengal		
Accreditation Standard	ISO/IEC 17025: 2005		
Discipline	Forensic Testing	Issue Date	20.11.2015
Certificate Number	T-1576	Valid Until	19.11.2017
Last Amended on	-	Page	5 of 17

S. No.	Product / Material of Test	Specific Test Performed	Test Method Specification against which tests are performed	Range of Testing / Limits of Detection
2.	Mobile phones, SIM cards, Memory cards & Smart phones (Blackberry, Android)	Data extraction of IMEI, Call logs, Contacts, Internet artifacts, E-mail GPRS settings, CDR	Working Procedure Manual in Computer Forensic CFSL/KOL/CF/WPM Issue No. 01; Date of Issue: 01.01.2007	Qualitative
3.	Secondary storage devices (CD, DVD, Pen Drive, Memory Card)	Data extraction and tracing origin of data to devices / users retrieval of detected / hidden data.	Working Procedure Manual in Computer Forensic CFSL/KOL/CF/WPM Issue No. 01; Date of Issue: 01.01.2007	Qualitative
4.	Devices with Internal Memory (Music Player, Digital Camera, Spycam, Medical Equipment)	Data extraction and tracing origin of data to devices / users retrieval of detected / hidden data	Working Procedure Manual in Computer Forensic CFSL/KOL/CF/WPM Issue No. 01; Date of Issue: 01.01.2007	Qualitative
5.	Data acquired from live memory and servers	Data extraction and tracing origin of data to devices / users retrieval of detected / hidden data	Working Procedure Manual in Computer Forensic CFSL/KOL/CF/WPM Issue No. 01; Date of Issue: 01.01.2007	Qualitative
6.	CCTV analysis	Data extraction / retrieval of detected / hidden data / timeline examination.	Working Procedure Manual in Computer Forensic CFSL/KOL/CF/WPM Issue No. 01; Date of Issue: 01.01.2007	Qualitative

Laboratory	Central Forensic Science Laboratory, 30, Gorachand Road, Kolkata, West Bengal		d, Kolkata,
Accreditation Standard	ISO/IEC 17025: 2005		
Discipline	Forensic Testing	Issue Date	20.11.2015
Certificate Number	T-1576	Valid Until	19.11.2017
Last Amended on	-	Page	6 of 17

S. No.	Product / Material of Test	Specific Test Performed	Test Method Specification against which tests are performed	Range of Testing / Limits of Detection
IV.	CHEMISTRY			
1.	Hand wash / Pocket wash /	Detection of Phenolpthalein		
	Swab from currency note in bribery cases	РН	CFSL/KOL/CHEM/WPM/02 Page. 11, Issue No. 01 Date of Issue: 01.01.2007	0 to 14
		Acid & Alkali	CFSL/KOL/CHEM/WPM/02 Page. 11; Issue No. 01 Date of Issue:01.01.2007	Qualitative
		Identification	CFSL/KOL/CHEM/WPM/02 Page. 11; Issue No. 01 Date of Issue:01.01.2007 (Thin Layer Chromatography, UV-Visible Spectrophotometry)	Qualitative Up to 1 μg/ml
		Detection of Carbonate ion	2 · · · · · · · · · · · · · · · · · · ·	
		Barium chloride	CFSL/KOL/CHEM/WPM/02 Page. 15; Issue No. 01; Date of Issue: 01.01.2007	1 mg/l
		Silver Nitrate test	CFSL/KOL/CHEM/WPM/02 Page. 15; Issue No. 01; Date of Issue:01.01.2007	1 mg/l
		Lime Water test	CFSL/KOL/CHEM/WPM/02 Page. 15; Issue No. 01; Date of Issue: 01.01.2007	Qualitative

Laboratory	Central Forensic Science Laboratory, 30, Gorachand Road, Kolkata, West Bengal		
Accreditation Standard	ISO/IEC 17025: 2005		
Discipline	Forensic Testing	Issue Date	20.11.2015
Certificate Number	T-1576	Valid Until	19.11.2017
Last Amended on	-	Page	7 of 17

S. No.	Product / Material of Test	Specific Test Performed	Test Method Specification against which tests are performed	Range of Testing / Limits of Detection
	Hand wash / Pocket wash / Swab from currency note in	Magnesium sulphate test	CFSL/KOL/CHEM/WPM/02 Page. 15; Issue No. 01; Date of Issue: 01.01.2007	1 mg/l
	bribery cases	Detection of Sodium ion		
		Flame test	CFSL/KOL/CHEM/WPM/02 Page. 15; Issue No. 01; Date of Issue: 01.01.2007	1 mg/l
		Zinc Uranyl Acetate Test	CFSL/KOL/CHEM/WPM/02 Page. 16; Issue No. 01; Date of Issue:01.01.2007	1 mg/l
V.	NARCOTICS			
1.	Plant Materials: Opium Morphine, Codeine	Marquis Test	CFSL/KOL/CHEM/WPM/01 Issue No. 01; Date of Issue: 01.01.2007	0.1 g/l
	Thebaine	Mecke's Test	CFSL/KOL/CHEM/WPM/01 Issue No. 01; Date of Issue: 01.01.2007	0.1 g/l
		Identification	CFSL/KOL/CHEM/WPM/01 Issue No. 01 Date of Issue: 01.01.2007 (Thin Layer Chromatography, Gas Chromatography –Mass Spectrometry)	0.1 g/l to 1.0 μg/l

Laboratory	Central Forensic Science Laboratory, 30, Gorachand Road, Kolkata, West Bengal		
Accreditation Standard	ISO/IEC 17025: 2005		
Discipline	Forensic Testing	Issue Date	20.11.2015
Certificate Number	T-1576	Valid Until	19.11.2017
Last Amended on	-	Page	8 of 17

S. No.	Product / Material of Test	Specific Test Performed	Test Method Specification against which tests are performed	Range of Testing / Limits of Detection
2.	Cannabis Ganja Charas, Hasish	Visualization for Morphological analysis (Macroscopic & Microscopic)	CFSL/KOL/CHEM/WPM/01 Issue No. 01 Date of Issue: 01.01.2007	Qualitative
		Color Test (Fast Blue B salt test)	CFSL/KOL/CHEM/WPM/01 Issue No. 01 Date of Issue: 01.01.2007	0.1g/l
		Identification	CFSL/KOL/CHEM/WPM/01 Issue No. 01 Date of Issue: 01.01.2007 (Thin Layer Chromatography/ Gas Chromatography –Mass Spectrometry)	0.1 g/l to 1.0 μg/l
3.	Drug of abuse: (Powder, tablets, capsules, street drugs samples <u>Benzodiazepines</u> (Nitrazepam,	Color test (Formaldehyde Sulphuric acid)	Clark's Isolation and Identification of Drugs, 2nd Edition, Page 134 CFSL/KOL/CHEM/WPM/02 Issue No. 01 Date of Issue: 01.01.2007	0.1g/l
	Alprazolam, Diazepam, Oxazepam, Clonazepam, Lorazepam	Identification	CFSL/KOL/CHEM/WPM/02 Issue No. 01 Date of Issue: 01.01.2007 (Gas Chromatography –Mass Spectrometry/Infra-red Spectroscopy)	1.0 μg/l

Laboratory	Central Forensic Science Laboratory, 30, Gorachand Road, Kolkata, West Bengal		
Accreditation Standard	ISO/IEC 17025: 2005		
Discipline	Forensic Testing	Issue Date	20.11.2015
Certificate Number	T-1576	Valid Until	19.11.2017
Last Amended on	-	Page	9 of 17

S. No.	Product / Material of Test	Specific Test Performed	Test Method Specification against which tests are performed	Range of Testing / Limits of Detection
4.	Barbiturates (Phenobarbital)	Color test (Dille-Koppanyi )	CFSL/KOL/CHEM/WPM/01 Issue No. 01 Date of Issue: 01.01.2007	0.1g/l
		Identification	CFSL/KOL/CHEM/WPM/01; Issue No. 01 Date of Issue: 01.01.2007 (Thin Layer Chromatography/ UV-Visible Spectrophotometry/ Infra-red Spectroscopy)	0.1g/l to1.0 μg/l
5.	Other Narcotic Drugs and Psychotropic substances	Identification	CFSL/KOL/CHEM/WPM/01 ; Issue No. 01 Date of Issue: 01.01.2007 (Thin Layer Chromatography/ Gas Chromatography –Mass Spectrometry/ UV-Spectroscopy Infra Red Spectroscopy)	0.1g/l to 1.0 μg/l
VI.	QUESTIONED DO	DCUMENTS	mila Red Specifoscopy)	
1.	Handwritings in originals & in reproduced copies	Physical Examination Microscopic Examination, Physical Matching, Use of different light sources in examination of alterations	Working Procedure Manual No. CFSL/KOL/DOCU/WPM Issue No. 02, Date of Issue. 16-08-2012 Chapter 6 to 8 & ASTM Standards E-2290-07A	Qualitative
2.	Alteration in Documents		Working Procedure Manual No. CFSL/KOL/DOCU/WPM Issue No. 02, Date of Issue. 16-08-2012 Chapter 9 & ASTM Standards E-2290-03	Qualitative

Laboratory	Central Forensic Science Laboratory, 30, Gorachand Road, Kolkata, West Bengal		
Accreditation Standard	ISO/IEC 17025: 2005		
Discipline	Forensic Testing	Issue Date	20.11.2015
Certificate Number	T-1576	Valid Until	19.11.2017
Last Amended on	-	Page	10 of 17

S. No.	Product / Material of Test	Specific Test Performed	Test Method Specification against which tests are performed	Range of Testing / Limits of Detection
3.	Printed Materials		Working Procedure Manual No. CFSL/KOL/DOCU/WPM Issue No. 02, Date of Issue. 16-08-2012 Chapter 10 & ASTM Standards E-2389-05	Qualitative
4.	Type written Materials		Working Procedure Manual No. CFSL/KOL/DOCU/WPM Issue No. 02, Date of Issue. 16-08-2012 Chapter 12 & ASTM Standards E-2494-07	Qualitative
5.	Stamp Impressions		Working Procedure Manual No. CFSL/KOL/DOCU/WPM Issue No. 02, Date of Issue. 16-08-2012 Chapter 16 to 17 & ASTM Standards E-2286-07, E2289-03	Qualitative
6.	Paper		Working Procedure Manual No. CFSL/KOL/DOCU/WPM Issue No. 02, Date of Issue – 16- 08-2012 Chapter 15 & ASTM Standards E-2325-05	Qualitative
7.	Ink		Working Procedure Manual No. CFSL/KOL/DOCU/WPM Issue No. 02, Date of Issue. 16-08-2012 Chapter 18 ASTM Standards E-1789-04	Qualitative
8.	Toner in Documents produced by DTP Process		Working Procedure Manual No. CFSL/KOL/DOCU/WPM Issue No. 02, Date of Issue. 16-08-2012 Chapter 13 to 14 & ASTM Standards E-2390-06, E2389-05	Qualitative

Laboratory	Central Forensic Science Laboratory, 30, Gorachand Road, Kolkata, West Bengal		
Accreditation Standard	ISO/IEC 17025: 2005		
Discipline	Forensic Testing	Issue Date	20.11.2015
Certificate Number	T-1576	Valid Until	19.11.2017
Last Amended on	-	Page	11 of 17

S. No.	Product / Material of Test	Specific Test Performed	Test Method Specification against which tests are performed	Range of Testing / Limits of Detection
VII.	EXPLOSIVES			
1.	Solid / liquid substances, Solid Powder, debris materials and post explosion residues for detection of Explosives , cations, anions and elements tested for			
	Nitrate ion	Issue No. 01; Date of Issue:		3.33mg/ml
	Sulphide ion	Sodium Nitroprusside Test	01.01.2007	3.33mg/ml
	Potassium ion	Sodium Cobalt nitrite Test		0.0333 mg/ml
	Ammonium ion	Nessler's Reagent Test		3.33 x 10-6 mg/ml
I	Arsenic ion	Reinsch's Test		0.033 mg/ml
	Sulphur (Pyridine extract)	Pyridine + NaOH Test		0.333 mg/ml

Aluminium Alizarine-S Test (NaOH extract)

0.333 mg/ml

Laboratory	Central Forensic Science Laboratory, 30, Gorachand Road, Kolkata, West Bengal		
Accreditation Standard	ISO/IEC 17025: 2005		
Discipline	Forensic Testing	Issue Date	20.11.2015
Certificate Number	T-1576	Valid Until	19.11.2017
Last Amended on	-	Page	12 of 17

S. No.	Product / Material of Test	Specific Test Performed	Test Method Specification against which tests are performed	Range of Testing / Limits of Detection
2.	Organic Explosive	s Components		
a.	TNT	Colour Test (Alcoholic KOH)	CFSL/KOL/EXP/WPM Issue No. 01; Date of Issue: 01.01.2007	Qualitative
		Identification	CFSL/KOL/EXP/WPM Issue No. 01; Date of Issue: 01.01.2007 (Thin Layer Chromatography/ HPLC/Gas Chromatography – Mass Spectrometry)	m/z ions:210,89,63, 30,134,193,180,149,76, 62.
b.	RDX	Color test (Thymol + H2SO4+Ethanol )	CFSL/KOL/EXP/WPM Issue No. 01; Date of Issue: 01.01.2007	Qualitative
		Identification	CFSL/KOL/EXP/WPM Issue No. 01; Date of Issue: 01.01.2007 (Thin Layer Chromatography/ HPLC/Gas Chromatography – Mass Spectrometry)	m/z ions: 46,30,120,42,75,44,128, 56,71,41.
c.	PETN	Color test (Di phenyl amine)	CFSL/KOL/EXP/WPM Issue No. 01; Date of Issue: 01.01.2007	Qualitative
		Identification	CFSL/KOL/EXP/WPM Issue No. 01; Date of Issue: 01.01.2007 (Thin Layer Chromatography/ HPLC/Gas Chromatography – Mass Spectrometry)	m/z ions: 76,46,57,55,56,60,43,47, 97,71.

Laboratory	Central Forensic Science Laboratory, 30, Gorachand Road, Kolkata, West Bengal		
Accreditation Standard	ISO/IEC 17025: 2005		
Discipline	Forensic Testing	Issue Date	20.11.2015
Certificate Number	T-1576	Valid Until	19.11.2017
Last Amended on	-	Page	13 of 17

S. No.	Product / Material of Test	Specific Test Performed	Test Method Specification against which tests are performed	Range of Testing / Limits of Detection
d.	TETRYL	Colour Test (Alcoholic KOH)	CFSL/KOL/EXP/WPM Issue No. 01; Date of Issue: 01.01.2007	Qualitative
		Identification	CFSL/KOL/EXP/WPM Issue No. 01; Date of Issue: 01.01.2007 (Thin Layer Chromatography/ HPLC/Gas Chromatography – Mass Spectrometry)	m/z ions: 30, 241, 75, 46, 91, 181, 77, 76, 90, 62
е.	NG	Colour Test (Alcoholic KOH)	CFSL/KOL/EXP/WPM Issue No. 01; Date of Issue: 01.01.2007	Qualitative
		Identification	CFSL/KOL/EXP/WPM Issue No. 01; Date of Issue: 01.01.2007 (Thin Layer Chromatography/ HPLC/Gas Chromatography – Mass Spectrometry)	m/z ions: 46, 76, 30, 58, 44, 43, 31, 117, 64, 89.
VIII.	PHYSICS			
1.	Tool marks (metal, fabrics,	Physical comparison	Working Procedure Manual Physics WPM – Physics-1	Qualitative
	wood)	Photography	Issue No. 01 Date of Issue: 01.01.2007	
		Microscopic examination		
		Test impression comparison		

Laboratory	Central Forensic Science Laboratory, 30, Gorachand Road, Kolkata, West Bengal		
Accreditation Standard	ISO/IEC 17025: 2005		
Discipline	Forensic Testing	Issue Date	20.11.2015
Certificate Number	T-1576	Valid Until	19.11.2017
Last Amended on	-	Page	14 of 17

S. No.	Product / Material of Test	Specific Test Performed	Test Method Specification against which tests are performed	Range of Testing / Limits of Detection
2.	Soil (Identification of	Physical comparison	Working Procedure Manual Physics WPM – Physics-1	Qualitative
	source correspondence	Photography	Issue No. 01 Date of Issue: 01.01.2007	
	comparison)	Microscopic examination Particle size distribution	2 01 1.5 0110112007	
		Density Gradient Analysis		
3.	Broken objects	Physical comparison	Working Procedure Manual Physics WPM – Physics-1	Qualitative
	(Broken parts of material like wood, metal, glass)	Photography	Issue No. 01 Date of Issue: 01.01.2007	
		Microscopic examination		
		Physical Matching		
4.	(Lottery Ticket, CounterfeitPhysics WPI PhotographyIssue No. 01	Physical comparison	Working Procedure Manual Physics WPM – Physics-1 Issue No. 01 Date of Issue: 01.01.2007	Qualitative
		Photography		
		Date of Issue. 01.01.2007		
		Fluorescence		
		Optical character analysis		
5.	Restoration of erased	Photography	Working Procedure Manual	Qualitative
	erased identification marks	Microscopic examination (Chemical Etching)	Physics WPM – Physics-1 Issue No. 01 Date of Issue: 01.01.2007	

Laboratory	Central Forensic Science Laboratory, 30, Gorachand Road, Kolkata, West Bengal		
Accreditation Standard	ISO/IEC 17025: 2005		
Discipline	Forensic Testing	Issue Date	20.11.2015
Certificate Number	T-1576	Valid Until	19.11.2017
Last Amended on	-	Page	15 of 17

Material of Test against which tests are Limits of Detection performed
--

## IX. TOXICOLOGY

1.	Viscera sample, blood, gastric lavage & poisons containing Methyl alcohol	Color test (Chromotropic acid)	Working Procedure Manual No. CFSL/KOL/TOX/WPM Issue No 01, Date of Issue: 01.01.2007	2 mg/ml
2.	Ethyl alcohol	Identification	Working Procedure Manual No. CFSL/KOL/TOX/WPM Issue No 01, Date of Issue: 01.01.2007 (Head space Gas chromatography-Mass Spectrometry)	2 µg/ml
		Iodoform test	Working Procedure Manual No. CFSL/KOL/TOX/WPM	2 mg/ml
		Sulphomolybdic acid test	Issue No 01, Date of Issue: 01.01.2007	2 mg/ml
		Identification	Working Procedure Manual No. CFSL/KOL/TOX/WPM Issue No 01, Date of Issue: 01.01.2007 (Head space Gas chromatography-Mass Spectrometry	2 µg/ml

Laboratory	Central Forensic Science Laboratory, 30, Gorachand Road, Kolkata, West Bengal		
Accreditation Standard	ISO/IEC 17025: 2005		
Discipline	Forensic Testing	Issue Date	20.11.2015
Certificate Number	T-1576	Valid Until	19.11.2017
Last Amended on	-	Page	16 of 17

S. No.	Product / Material of Test	Specific Test Performed	Test Method Specification against which tests are performed	Range of Testing / Limits of Detection
3.	Phosphide	Color test (Mercuric Chloride)	Working Procedure Manual No. CFSL/KOL/TOX/WPM Issue No 01, Date of Issue: 01.01.2007	2 µg/ml
4.	Cyanide	Color test (Ferric thiocyanate )	Working Procedure Manual No. CFSL/KOL/TOX/WPM Issue No 01, Date of Issue: 01.01.2007	2 mg/ml
5.	Viscera sample, blood, gastric lavage having Arsenic, Mercury, Bismuth, Copper, Antimony, Lead, Aluminium, Zinc,	Reinsch Test, Potassium Chromate; Sodium Stannite; Ammonium hydroxide Test, Rhodamine-B; Dithizone Test, Dil. Nitric acid test, Aluminone test, Alizarin Test & Potassium ferrocyanide test	Working Procedure Manual No. CFSL/KOL/TOX/WPM Issue No 01, Date of Issue: 01.01.2007)	2 mg/ml
	Nitrate, Bromide, Iodide, Sulphate, Chromate, Flouride & Chloride	Copper Sulphate Test, Sulphanalic acid test, Silver Nitrate test, Lead acetate test & Barium Chloride test.	Working Procedure Manual No. CFSL/KOL/TOX/WPM Issue No 01, Date of Issue: 01.01.2007	2 mg/ml
		Identification	Working Procedure Manual No. CFSL/KOL/TOX/WPM Issue No 01, Date of Issue: 01.01.2007 (ICP-Mass spectrometry and Atomic Absorption Spectroscopy)	2 μg/ml

Laboratory	Central Forensic Science Laboratory, 30, Gorachand Road, Kolkata, West Bengal		
Accreditation Standard	ISO/IEC 17025: 2005		
Discipline	Forensic Testing	Issue Date	20.11.2015
Certificate Number	T-1576	Valid Until	19.11.2017
Last Amended on	-	Page	17 of 17

S. No.	Product / Material of Test	Specific Test Performed	Test Method Specification against which tests are performed	Range of Testing / Limits of Detection
6.	Viscera sample, blood, gastric lavage, containing Insecticides. Malathion,	Color test (Fast Blue B salt & Tollens reagent)	Working Procedure Manual No. CFSL/KOL/TOX/WPM Issue No 01, Date of Issue: 01.01.2007	2 mg/ml
	parathion, Endosulfan, Iso-propoxur Allethrin & Cypermethrin	Identification	Working Procedure Manual No. CFSL/KOL/TOX/WPM Issue No 01, Date of Issue: 01.01.2007 (Thin layer chromatography/ Gas chromatographic -mass spectrometry)	0.1 g/l to 1.0 μg/l
7.	Viscera sample, blood, gastric lavage containing alkaloids , Benzodiazepines,	Color test (Dragon droffs & Forrest reagent.)	Working Procedure Manual No. CFSL/KOL/TOX/WPM Issue No 01, Date of Issue: 01.01.2007	2 mg/ml
	Barbiturates, Paracetamol	Identification	Working Procedure Manual No. CFSL/KOL/TOX/WPM Issue No 01, Date of Issue: 01.01.2007 (Thin layer Chromatography/ Gas chromatography-mass spectrometry)	0.1 g/l to 1.0 μg/l
8.	Viscera sample, blood, gastric lavage & poisons containing glass powder, diamond powder, hair	Microscopic examination & Flame test	Working Procedure Manual No. CFSL/KOL/TOX/WPM Issue No 01, Date of Issue: 01.01.2007	Qualitative