

Laboratory	Devansh Testing & Research Laboratory Pvt. Ltd., 94, Shiv Ganga Industrial Estate, Lakeshari, Bhagwanpur, Roorkee, Haridwar, Uttarakhand		
Accreditation Standard	ISO/IEC 17025: 2005		
Discipline	Biological Testing	Issue Date	24.11.2013
Certificate Number	T-2089	Valid Until	23.11.2015
Last Amended on	-	Page	1 of 4

S.No.	Product / Material of Test	Specific Test Performed	Test Method Specification against which tests are performed	Range of Testing / Limits of Detection
--------------	-----------------------------------	--------------------------------	--	---

I. DRUGS & PHARMACEUTICALS

A. Raw Materials & Finish Products

1.	Microbial Contamination of Non-Sterile Products (Maize starch, Purified water, Tablets etc.)	Total Aerobic Microbial Count	IP 2010 Page No.46 Method (2.2.9)	1/10 cfu to 10000 cfu / ml or gm. or as specified
		Total yeast & Molds	IP 2010 Page No.46 Method (2.2.9)	(Present / Absent)/ gm or as specified
		Escherichia coli	IP 2010 Page No.46 Method (2.2.9)	(Present / Absent) / gm or as specified
		Salmonella Sp.	IP 2010 Page No.46 Method (2.2.9)	(Present / Absent) / gm or as specified
		Staphylococcus aureus	IP 2010 Page No.46 Method (2.2.9)	(Present / Absent) / gm or as specified
		Pseudomonas aeruginosa	IP 2010 Page No.46 Method (2.2.9)	(Present / Absent) / gm or as specified
2.	Gentamycin Sulphate I.P.	Assay	IP 2010 Page No. 49 Method A (2.2.10)	± 20 % of specified limits
3.	Amikacin Sulphate I.P.	Assay	IP 2010 Page No. 49 Method B (2.2.10)	± 20 % of specified limits

II. FOOD & AGRICULTURAL PRODUCTS

1.	Milk Powder IS 1165: 2002; RA 2009	Total Plate Count	IS 5402:2012	10 to 100000 cfu/gm
		Coliform	IS 5401 (P-1):2012	Absent/Present per 0.1gm
		S.aureus	IS 5887 (P-2): 1976 RA2009	Absent/Present per 0.1gm

Laboratory	Devansh Testing & Research Laboratory Pvt. Ltd., 94, Shiv Ganga Industrial Estate, Lakeshari, Bhagwanpur, Roorkee, Haridwar, Uttarakhand		
Accreditation Standard	ISO/IEC 17025: 2005		
Discipline	Biological Testing	Issue Date	24.11.2013
Certificate Number	T-2089	Valid Until	23.11.2015
Last Amended on	-	Page	2 of 4

S.No.	Product / Material of Test	Specific Test Performed	Test Method Specification against which tests are performed	Range of Testing / Limits of Detection
		Salmonella	IS 5887 (P-3):1999; RA2009	Absent/Present per 25gm
		Shigella	IS 5887 (P-7):1999; RA2009	Absent/Present per 25gm
2.	Cream (Pasteurized Cream)	Total Plate Count	IS 3509:1966; RA2013	10 to 100000 cfu/gm
		Coliform Count	IS 3509:1966; RA2013	10 to 1000 cfu/gm
3.	Butter	Coliform Count	IS 3507:1966; RA2013	10 to 1000 cfu/ml
		Yeast & Mould Count	IS 3507:1966; RA2013	10 to 1000 cfu/ml
4.	Malted Milk Food (IS 1806: 1975)	Total Plate Count	IS 5402:2011	10 to 100000 cfu/gm
		Coliform Count	IS 5401 (P-1):2011	10 to 1000 cfu/gm
5.	Dahi IS 9617: 1980	Yeast & Mould Count	IS 5403:1999; RA2009	10 to 1000 cfu/gm
		Coliform	IS 5401 (P-1):2009	10 to 1000 cfu/ gm
6.	Milk (Pasteurized Milk) IS 1479 (P-3) 1977 RA2009	Total Plate Count	IS 5402:2012	1 to 100000 cfu/ml
		Coliform	IS 5401 (P-1): 2012	Absent/Present per /ml
		S.aureus	IS 5887 (P-2): 1976; RA2009	Absent/Present per ml
7.	Condensed Milk IS 1166: 1986; RA1997	Total Plate Count	IS 5402: 2002	10 to 10000 cfu/gm
		Coliform	IS 5401 (P-1):2002	Absent/Present per 0.1 gm
		Yeast & Mould Count	IS 5403:1999; RA2005	10 to 100 cfu/ gm
		S.aureus	IS 5887 (P-2):1976; RA2005	Absent/Present per gm

Laboratory	Devansh Testing & Research Laboratory Pvt. Ltd., 94, Shiv Ganga Industrial Estate, Lakeshari, Bhagwanpur, Roorkee, Haridwar, Uttarakhand		
Accreditation Standard	ISO/IEC 17025: 2005		
Discipline	Biological Testing	Issue Date	24.11.2013
Certificate Number	T-2089	Valid Until	23.11.2015
Last Amended on	-	Page	3 of 4

S.No.	Product / Material of Test	Specific Test Performed	Test Method Specification against which tests are performed	Range of Testing / Limits of Detection
8.	Milk Cereal Based Food IS 1656: 2007; RA 2009	Total Plate Count	IS 5402:2002	10 to 10000 cfu/gm
		Coliform	IS 5401 (P-1): 2002	Absent/Present per 0.1 gm
		E.Coli	IS 5887 (P-1): 1976; RA2005	Absent/Present per 0.1 gm
		S.aureus	IS 5887 (P-2): 1976; RA2005	Absent/Present per 0.1 gm
		Salmonella	IS 5887 (P-3): 1976; RA2005	Absent/Present per 25 gm
		Shigella	IS 5887 (P-7): 1976; RA2005	Absent/Present per 25 gm
		Yeast & Mould Count	IS 5403 : 1999; RA2005	Absent/Present per 0.1 gm
9.	Chocolates IS 1163: 1992; RA 2004	S.aureus	IS 5887 (P-2):1976; RA2005	Absent/Present per 10 gm
		Salmonella	IS 5887 (P-3):1999; RA2005	Absent/Present per 25 gm
		E.Coli	IS 5887 (P-1):1999; Reaff:2005	Absent/Present per 10 gm
		Yeast & Mould Count	IS 5403 : 1999; RA2005	10 to 1000 cfu/gm
10.	Corn Flour IS 1005: 1992; RA2009	E.Coli	IS 5887 (P-1): 1976; RA2005	Absent/Present per 25 gm
		Salmonella	IS 5887 (P-3) : 1999; A2005	Absent/Present per 25 gm
11.	Non-Alcoholic Beverages IS 13019: 1991	Total Plate Count	IS 5402:2002	1 to 1000 cfu/ml
		Coliform	4831:2006 / (IS 5401 (P-2) 2002)	< 0.3 to >110 MPN/100ml
		Yeast & Mould Count	IS 5403:1999; RA2005	1 to 500 cfu/ml

Laboratory **Devansh Testing & Research Laboratory Pvt. Ltd., 94, Shiv Ganga Industrial Estate, Lakeshari, Bhagwanpur, Roorkee, Haridwar, Uttarakhand**

Accreditation Standard **ISO/IEC 17025: 2005**

Discipline **Biological Testing** **Issue Date** **24.11.2013**

Certificate Number **T-2089** **Valid Until** **23.11.2015**

Last Amended on **-** **Page** **4 of 4**

S.No.	Product / Material of Test	Specific Test Performed	Test Method Specification against which tests are performed	Range of Testing / Limits of Detection
12.	Carbonated Beverages IS 2346: 1992; RA 2009	Total Plate Count	IS 5402:2002	1 to 1000 cfu/ml
		Coliform	ISO 4831:2006/ IS 5401 (P-2) 2002	Absent/Present per 100 ml
		Yeast & Mould Count	IS 5403:1999; RA2005	1 to 500 cfu/ml
-X-X-X-X-X-X-X-X-X-X-X-X-				